



Behind Every Smile is a Story

- QUALITY - ЯКІСТЬ
- PRECISION - ТОЧНІСТЬ
- ACCURACY - ПРАВИЛЬНІСТЬ
- DURABILITY - ЗНОСОСТІЙКІСТЬ



# ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтичні
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**MD-1800**  
Height Gauge 0.18  
Aluminum



**MD-1801**  
Height Gauge 0.22  
Aluminum



**MD-1802**  
Height Gauge 0.18  
Stainless Steel



**MD-1803**  
Height Gauge 0.22  
Stainless Steel



**MD-1804**  
Moveable  
Height Gauge



**MD-1805**  
Endo Finger  
Ruler Span  
Measure Scale



**MD-1806**  
Endo Finger  
Ruler Span  
Measure Scale



**MD-1807**  
Endo Finger  
Ruler Span  
Measure Scale



**MD-1808**  
Positive  
Postioner



**MD-1809**  
Strip Saw



**MD-1810**  
Bone Bracket Gauge



**MD-1811**  
Endo Finger Ruler  
Span Measure Scale



**MD-1812**  
Endo Finger Ruler  
Span Measure Scale

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**MD-1813**  
Measuring Loop  
angled, length  
175 mm,  
measuring range  
0 - 40 mm



**MD-1814**  
Measuring Loop  
straight length  
175 mm,  
measuring range  
0 - 40 mm



**MD-1815**  
Measuring loop, Castroviejo  
angled, length 175 mm,  
measuring range 0 - 40 mm,  
with window



**MD-1816**  
Measuring Loop  
straight length 90 mm,  
measuring range  
0 - 40 mm,

**MD-1817**  
Castroviejo Calipers  
20mm  
90mm



**MD-1818**  
Measuring Loop  
angled, length 90 mm,  
measuring range  
0 - 40 mm,

**MD-1819**  
Castroviejo Calipers  
20mm, 45°  
90mm



**MD-1820**  
Castroviejo Calipers  
20mm  
90mm



**MD-1821**  
Fig. 1  
Iwanson  
spring caliper  
For metal and  
porcelain.



**MD-1822**  
Fig. 2  
Iwanson  
spring caliper  
For wax.

# ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**MD-1823**  
Boley Gauge



**MD-1824**  
Vernier Gauge



**MD-1825**  
Aluminium Caliper  
up to 200mm



**MD-1826**  
Boley Gauge



**MD-1827**  
Delux Boley  
Gauge



**MD-1828**  
Beerendonk  
130mm  
Mini Gauge



**MD-1829**  
Compasso  
De Willis



**MD-1830**  
Dentaurum  
Ruler



**MD-1831**  
Ruler Wedge



**MD-1832**  
Ruler

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтичні
- Instruments Orthodontiques
- أدوات تقويم الأسنان

**Bracket Holder Tweezers**



**MD-1833**  
Deluxe Bracket  
Holder Tweezers



**MD-1834**  
Deluxe Bracket  
Holder Tweezers



**MD-1835**  
Direct Bond Bracket  
Holder Tweezers



**MD-1836**  
Bracket Tweezers  
Slot Aligner



**MD-1837**  
Direct Bond  
Bracket Holder  
Tweezers, Small Tip



**MD-1838**  
Exact Bracket  
Holder Tweezers



**MD-1839**  
Bracket Placing  
Tweezers



**MD-1840**  
Easy Access  
Tweezers



**MD-1841**  
Bracket Placer  
Tweezers  
(Extra Fine Tip)



**MD-1842**  
Deluxe Bracket  
Holder Tweezers  
135mm



**MD-1843**  
Easy Eccess  
Tweezers



**MD-1844**  
Slipper  
Tweezers



**MD-1845**  
Direct Bonding  
Bracket Holders  
135mm



**MD-1846**  
Direct Bonding  
Bracket Holders  
136mm



**MD-1847**  
Direct Bonding  
Bracket Holders  
135mm



**MD-1848**  
Deluxe Bracket  
Holder



**MD-1849**  
Direct Bonding  
Bracket Holders  
140mm

# ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- ادوات تقويم الأسنان



**MD-1850**  
Band seater  
auto cleavable



**MD-1851**  
Band seater  
auto cleavable  
Triangle Shape



**MD-1852**  
Band seater  
auto cleavable  
Square Shape



**MD-1853**  
Band seater  
auto cleavable  
Round Shape



**MD-1854**  
Mershon Band  
Pusher  
140mm



**MD-1855**  
Band seating &  
removing  
instruments



**MD-1856**  
Ortho Light double  
ended  
ligature director  
scaler



**MD-1857**  
Ortho Light Band  
Pusher Scaler



**MD-1858**  
Band Seating &  
Removing  
Silicon Grip



**MD-1859**  
Ligature Instrument  
Silicon Grip



**MD-1860**  
Elastic Instrument  
Silicon Grip



**MD-1861**  
Distal Bender  
Silicon Grip



**MD-1862**  
Distal Bender  
Silicon Grip



**MD-1863**  
Distal Bender  
Silicon Grip

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**MD-1864**  
Ortho light scaler double ended



**MD-1865**  
Ortho light scaler double ended



**MD-1866**  
Ortho light scaler double ended



**MD-1867**  
Twiel-on elastic module instrument



**MD-1868**  
Twiel-on elastic module instrument



**MD-1869**  
Schure Band Seater



**MD-1870**  
Ortho light scaler



**MD-1871**  
Ortho light scaler



**MD-1872**  
Ortho light scaler heavy



**MD-1873**  
Ortho instrument



**MD-1874**  
Ortho instrument

# ORTHODONTIC INSTRUMENTS

- Инструменты ортодонтические
- Instrumentos Ortodonticos
- Orthodontische Instrumente
- Инструменти ортодонтични
- Instruments Orthodontiques
- أدوات تقويم الأسنان



**MD-1875**  
Ortho Light double ended  
ligature director



**MD-1876**  
Ortho Light double ended ligature scaler director tivate coil spring  
are remove modules with one instruments



**MD-1877**  
Ortho band pusher scaler



**MD-1878**  
Ortho ligature director single ended



**MD-1879**  
Ortho Light double ended ligature  
director scaler



**MD-1880**  
Ortho Light double ended ligature scaler/director  
tivate coil spring are remove modules with one instruments



**MD-1881**  
Ortho Light double ended ligature director



**MD-1882**  
Ortholite Ligature Director Scaler



**MD-1883**  
Ortho light ligature director double ended



**MD-1884**  
Double ended cement spatula



**MD-1885**  
Ortho light scaler single ended



- Щипцы для ортодонтии и протезирования
- Pincès pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستتیس

**BOX JOINT PLIERS**
**MD-1886**

Universal Dolphin Style  
Wire Bending & cutting pliers


**MD-1887 Band Removing Plier**

The curved blade of the band removing plier can strip easily and gently remove directly to the tooth. Length 130mm


**MD-1888**

Tweed Loop forming Pliers


**MD-1889**

**How Plier curved**  
Universally usable tool for the multiband technology for wires up to .028". Length 130mm


**MD-1890 How Plier**

Universally usable tool for the multiband technology for wires up to .028". Length 135mm


**MD-1891**

Nance Loop forming Pliers


**MD-1892**

**Bracket Removing Plier curved**  
Bracket Removing Plier curved for easy removal of brackets in the posterior region. Length 130mm


**MD-1893**

Young Wire Bending & Loop Forming Pliers


**MD-1894**

Wire Bending & Loop Forming Pliers


**MD-1895**


Adams, small wire bending and loop forming Pliers


**MD-1896**

**Bracket Removing Plier Round Box**  
Bracket removing plier for comfortable and quick removal of the brackets. Length 135mm


**MD-1897**


**Weingart Plier medium**  
Weingart plier that are ideal due to their curved tip for setting sheets and is capable of bending the ends. Length 145mm




- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستيتيس

## BOX JOINT PLIERS

**MD-1898 Adams Plier Round Box**  
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



**MD-1904 Adams Plier**  
Adams pliers for bending wire arches and bends right angles for wires up to .022", .025 tip to beak. Length 125mm



**MD-1899 Falcon Beak Contouring Pliers.**



**MD-1905 Tweed Loop Formin Pliers**



**MD-1900**  
wires up to 0.020" / 0.7mm.



**MD-1906 Bird Beak Plier Standard**  
Bird Beak Wire bending plier without notches, which is suitable as a universal utility tool. Length 125mm



**MD-1901 Johnson Contouring Pliers**



**MD-1907 Mini Ligature cutters**



**MD-1902 Aderer Large wire bending and loop forming pliers**



**MD-1908 Small Fluch-cut distal and cutting Pliers** with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



**MD-1903 Tweed / angle rectangular arch forming & Torquing pliers**



**MD-1909 Large Fluch-cut distal and cutting Pliers** with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips



- Щипцы для ортодонтии и протезирования
- Pincettes pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمشة لتقويم الأسنان بروستتیس

**BOX JOINT PLIERS**
**MD-1910**

**Tweed / Angel**  
 1.78mm (.070") wide tips arch  
 forming and torquing pliers


**MD-1916**

**Contouring Plier**


**MD-1911**

**Tweed / Angel**  
 1.78mm (.070") wide tips arch  
 forming and torquing pliers


**MD-1917**

**Bands & Crowns Plier**  
 in order to skirt


**MD-1912**

**Lingual / Palatal arch**  
 forming pliers


**MD-1918**

**Ligature Forming Pliers**


**MD-1913**

**Aderer**  
 Facebow loop forming Pliers


**MD-1919**

for double back and triple back  
 bands in 0.30" / 0.76mm & 0.36"  
 / 0.91mm lingual arch wires


**MD-1914**

**Waldsavhs** wire bending  
 and cutting pliers


**MD-1920**

**Bird Beak** Wire Bending  
 & cutting pliers


**MD-1915 Aderer Plier**

For the manufacturing of appliances and  
 activation of wires for wires up to .030" .fine  
 tip and up to .036" mean tip.  
 Length 130mm


**MD-1921**

**Weingart Plier medium**  
 Weingart plier that are ideal due to  
 their curved tip for setting sheets and  
 is capable of bending the ends.  
 Length 145mm



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستيثيس

## L-KEY JOINT PLIERS



### MD-1922

Distal end cutter with safety hold, mini head, long handles. fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



### MD-1923

Flush-cut distal end cutter with safety hold fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" to .025")



### MD-1924

Soft wire and ligature cutter, Mini fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



### MD-1925

Soft wire and ligature cutter, standard fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



### MD-1926

Soft wire and ligature cutter, 15° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



### MD-1927

Soft wire and ligature cutter, 45° angled fitted with un-breakable, corrosion resistant and durable stainless steel tips

- ☞ .25 mm to .38 mm (.010" to .015")



### MD-1928

Hard wire cutter, straight Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### MD-1929

Hard wire cutter, 15° angled Fitted with un-breakable, corrosion resistant and durable stainless steel tips

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### MD-1930

Hard wire cutter, 15° angled Non-autoclavable

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كاشاة لتقويم الأسنان بروستتیس

## L-KEY JOINT PLIERS



**MD-1931**  
Band removing pliers



**MD-1936**  
Weingart-mini  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")




**MD-1932**  
Bracket removing pliers angled



**MD-1937**  
Weingart  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



**MD-1933**  
Bracket removing pliers straight



**MD-1938**  
How Straight  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")




**MD-1934**  
Bracket removing pliers with plastic tip



**MD-1939**  
Angled  
Utility Pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



**MD-1935**  
Adhesive removing pliers



**MD-1940**  
Wire bending and loop forming pliers  
Angled (Bird beak)

- 0.30 mm to 0.50 mm (.012" to .020")
- Max 0.56 mm x 0.71 mm (.022" x .028")



# PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستيتيس

## L-KEY JOINT PLIERS



### MD-1941

Multi-purpose pliers, bird beak, De la rosa and cutter in one

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.56 mm (.018" x .028")



### MD-1946

Stop pliers V-Bebd - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.43 mm x 0.56 mm (.017" x .028")



### MD-1942

Cinch Back / tie back Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### MD-1947

Lingual / palatal arch forming plies

- 0.76 mm (.030") / 0.91 mm (.036")



### MD-1943

Hammer head Distal end bending pliers

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### MD-1948

Torquing pliers set



### MD-1949

Tweed / angel Rectangular arch forming and torquing pliers

- 0.30 mm to 0.71 mm (.012" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### MD-1944

Nance Loop forming pliers

- 0.36 mm to 0.71 mm (.014" to .028")
- Max. 0.56 mm x 0.71 mm (.022" x .028")



### MD-1950

Surgical ball hook crimping pliers, angled



### MD-1945

Detailing step pliers / Bayonet bending pliers - 1mm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")




### MD-1951

Surgical ball hook crimping pliers, straight




- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستيتيس


## L-KEY JOINT PLIERS




**MD-1952**  
Aderer Clamp three jaw wire bending plier for contouring and bending wire up to .30"/.76mm.




**MD-1958**  
Falcon beak contouring plier.





**MD-1953**  
Bird beak plier for bending round, square and rectangle wire up to .30"/.76mm.



**MD-1959**  
For placement and removal of arch wires, pins and other auxiliaries.




**MD-1954**  
Bird beak plier for light wire bending round, square and rectangular up to .016"/0.41mm.



**MD-1960**  
For placement and removal of arch wires, pins and other auxiliaries in angle design



**MD-1955**  
Double Rounded jaw bending pliers, hooks and wire up to .020"/0.51mm.





**MD-1961**  
contouring Plier.




**MD-1962**  
Bending plier. In order to skirt molares bands. With grooved. Recommended wire. 016", 018" and 022" 0.41mm, 0.46mm and 0.56mm



**MD-1956**  
Adams pliers for hooks and bending wire up to .030"/0.76mm.



**MD-1963**  
Loop forming in four length 3,4,5 & 6 mm wire capacity 0.22", 56mm.



**MD-1957**  
Johnson pliers for crown band contouring



**MD-1964**  
Arch bending plier design for convenient handling of square or rectangular wire, .020" / 0.51mm.



- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستیتیس


## LAP & L-KEY JOINT PLIERS



**MD-1965**  
In order to skirt bands and crowns.



**MD-1966**  
Wires up to .020" / 0.7mm.



**MD-1967**  
Round and concave plier for loop forming and banding wire up to .20" / .051mm.




**MD-1968**  
Bending for round square or rectangular wire up to 0.30"/0.76mm



**MD-1969**  
for double back and triple back bands in 0.30" / .76mm & 0.36" / .91mm lingual arch wires.



**MD-1970**  
Detailing step plier double-sided offset beak permits bayonet bends. Creates right hand or left hand step-up or stepdown at 1/2, 3/4, or 1mm increments.



**MD-1971**  
"V" bend stop special design makes 1mm "v" bends to shorten archwires or provide a positive stop. Excellent for placing stops in titanium wires. For round or rectangular wires up to .025"



**MD-1972**  
Angle how plier efficient utility plier with serrated beaks for a firm grip. 45 angled tips to access hard to reach areas.



**MD-1973**  
Weingardt plier serrated tips hold wires firmly at a convenient working angle. Tapered beaks fit easily between brackets and hard to reach areas. Our most popular utility plier.



**MD-1974**  
Pin & Ligature Cutters TC.



**MD-1975**  
Mini Pin & Ligature Cutter TC.



**MD-1976**  
Standard Ligature Cutter TC.



**MD-1977**  
V Bend Stop Plier TC.



**MD-1978**  
Angled Direct Bond Bracket Remover TC.



**MD-1979**  
Cutter Str.

Distal Cutter 14cm  
**MD-1980** 2.5 mm  
**MD-1981** 3.0 mm  
**MD-1982** 3.0 mm



**MD-1983**  
Arrow Clasp Forming Pliers 12.0cm

**MD-1984**  
Schwarz 12.5cm





**MD-1985**  
 Distal end Cutter, Safety Hold Plier  
 (Ideal for Cutting Wire) TC.



**MD-1986**  
 Hook Crimping Plier Angled  
 Straight TC.



**MD-1987**  
 3 Prong Clasp Adjustable  
 Plier TC.



**MD-1988**  
 Utility Plier Weingart TC.



**MD-1989**  
 Direct Bracket Remover  
 TC.



**MD-1990**  
 Straight Direct Bond Bracket  
 Remover TC.



**MD-1991**  
 Pin Holding Plier TC.



**MD-1992**  
 Hard wire cutters



**MD-1993**  
 Elastic Separating Plier  
 Round Handle



**MD-1994**  
 Elastic  
 Separating Plier




**MD-1995 Str.**  
**MD-1996 Cvd.**  
 Smaha  
 Elastomeric  
 ligature  
 placing  
 forceps  
 140 mm



**MD-1997**  
 Elastomeric  
 ligature  
 placing  
 forceps  
 with  
 groove  
 Falcon 140 mm



**MD-1998**  
 Elastomeric  
 ligature  
 placing  
 forceps  
 with  
 groove  
 Falcon 140 mm




- Щипцы для ортодонтии и протезирования
- Pincas pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروستیتیس

## CUTTER, UTILITY PLIERS & WIRE FORMING PLIERS




**MD-1999**  
Weingart Pliers



**MD-2004**  
Adhesive  
Removing Pliers




**MD-2000**  
Slim Weingart  
Pliers



**MD-2005**  
Posterior Band Removing  
Pliers - Long




**MD-2001**  
How Pliers  
Straight




**MD-2006**  
Posterior Band Removing  
Pliers - Short



**MD-2002**  
How Pliers Offset




**MD-2007**  
Bird Beak Pliers



**MD-2003**  
Band seating Pliers



Replacement tip and  
hex key included.



**MD-2008**  
Omega Loop  
Forming Pliers

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонтії та протезування
- Alicates para ortodoncia y prótesis
- كماشة لتقويم الأسنان بروستیتیس

## CUTTERS



**MD-2009**  
 Hard  
 Wire Cutter 15




**MD-2013**  
 Mini Pin &  
 Ligature Cutter




**MD-2010**  
 Distal End Cutter  
 Flush Cut




**MD-2014**  
 Hard Wire Cutter  
 Straight



**MD-2011**  
 Universal Cut & Hold  
 Distal End Cutter



**MD-2015**  
 Micro-Mini Pin &  
 Ligature Cutter



**MD-2012**  
 Pin & Ligature  
 Cutter




**MD-2016**  
 Cutter

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
- Pincers pour l'orthodontie et la prothèse
- Zangen für Kieferorthopädie und Prothetik
- Щипці для ортодонції та протезування
- Alicates para ortodoncia y prótesis
- كمانشة لتقويم الأسنان بروستيثيس



**MD-2017**  
Three Jaw Pliers




**MD-2021**  
Jarabak Pliers




**MD-2018**  
Light Wore Bird Beak  
Pliers with Groove




**MD-2022**  
Hollow Chop  
Contouring Pliers




**MD-2019**  
Tweed loop  
Forming Pliers



**MD-2023**  
Lingual Arch  
Forming Pliers




**MD-2020**  
Arch Bending Pliers



**MD-2024**  
Rectangular Arch  
Bending Pliers

## ORTHODONTICS PLIERS STANDS

- Ортодонтия Щипцы Подставки
- Supports de pinces d'orthodontie
- Kieferorthopädie Zangenständer
- Ортодонтия Щипці Підставки
- KFO Zangen Ständer
- تقويم الأسنان كماشة تقف



**MD-2025**  
Plier Stand  
Stainless Steel.



**MD-2026**  
Plier Stand  
Stainless Steel.



**MD-2027**  
Pliers Stand  
Stainless Rod  
(Aluminum Anodized)  
Stainless Steel



**MD-2028**  
Plier Stand  
Stainless Steel.



**MD-2029**  
Plier Stand  
Stainless Steel.



**MD-2030**  
Plier Stand  
Plastic